

*Exceptional service in the national interest*



# Analysis of Hyperspectral Images of Chlamy on Lidke Line Scanner

Imaging: Aysha and Keith

Analysis: Jeri

Data taken on 7/9/14, 7/10/14,  
7/16/14, 7/17/14



Sandia National Laboratories is a multi-program laboratory managed and operated by Sandia Corporation, a wholly owned subsidiary of Lockheed Martin Corporation, for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-AC04-94AL85000. SAND NO. 2011-XXXXP

# Analysis Details

- \*.hsi files were converted to s3d files using batch script that reads in data with Lidke code, then calls Creates3d.m. Expansion factor was utilized for improved precision.
- \*.s3d files were opened individually in ShowMe3D.exe and checked for integrity/quality.
- Sets of \*.s3d files that were to be compared were combined into a composite data set using <Preprocess> function in ImageMCR.
- MCR was run on resulting composite files, excluding any individual traces known to contain saturated pixels (more on this later – Slide 4).


# Observation 1

- The baseline does not appear to be a linear offset as we are used to. Instead it seems to slope upward toward the red wavelengths.  
Note: the data is often quite a bit below zero.
  - This could be adequately accounted for if we had some dark current images of the same size, but these were not present.
  - Instead a linear offset was modelled and allowed to be negative to improve the results.

**\*\* Future Recommendation:** Collect a dark image (no light on the detector at the beginning, middle and end of the each day. Include those in with the data files.

# Observation 2

- Many of the files exhibited the classic “rolled-over peak) indicative of detector nonlinearity. This was verified by opening raw \*.hsi files in hyperview (keith’s software).
  - This leads to strong broadening effect on data, and in most files is too severe to perform MCR analysis.



**\*\* Future Recommendation:** Choose a laser power or integration time to ensure no signal is greater than the linear range of detector. In this case I am guessing about 1500 raw cts, but Keith would know better.


# Observation 3

- The signal intensity varied greatly from cell to cell even within similar conditions
  - Occasionally this was due to the cell being dead
  - Sometimes it looked like cell wasn't well focused?
  - Was the laser power or integration time adjusted from image to image? If so this is not recommended.

**\*\* Future Recommendation:** Don't vary power or integration during an experiment. Improvements in focusing were definitely seen as you became more experienced!

# Encapsulated Chlamy


(Nadia bicarb experiment)




This 2<sup>nd</sup> derivative like shape represents broadening and is due to nonlinear effect of bright pixels.

Note: Spectral residuals are still quite structured indicating additional shifts and broadenings that are unaccounted for, but represent only a small portion of the variance with this model, therefore the model is deemed sufficient.

# Component 2 (CHL-685)



# Component 3 (LHCII?)






# Encapsulated Chlamy


Nadia – BiCarb Thick



Nadia – BiCarb Thin




Nadia – Media Thick



\*signal for this row was much weaker than others?

Nadia – Media Thin



Color scales are adjusted independently for each image.

# Summary Encapsulated Chlamy

(Nadia bicarb experiment)

- Not enough consistent data to determine if there are any trends.
  - There's some segmentation of the red and green signal, however it doesn't appear to trend with condition.
  - There could be an overall change in intensity, but it is not certain whether this is due to condition or other factors like focusing.

**\*\* Future Recommendation:**  
Continue collecting data. We typically like 30-50 cells at each condition for robust statistics.