

Exceptional service in the national interest



Analysis of Fixed Retina Samples

Sample preparation: Carolina Nitta

Imaging: Mike Sinclair (7-18-14)

Analysis: Jeri Timlin (7-27-14)





Sandia National Laboratories is a multi-program laboratory managed and operated by Sandia Corporation, a wholly owned subsidiary of Lockheed Martin Corporation, for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-AC04-94AL85000. SAND NO. 2011-XXXXP

Imaging/Analyiss Details

- 488 nm excitation, OD 1 used for GFP expressing samples, OD 0 used for non-GFP samples
- 20x objective, 100 x 100 μm field of view
- Spatial resolution: approximately 500nm in X&Y and 1.0 μm in Z
- Analysis: MCR analysis (ImageMCR) was run on a composite set of all the images taken (compressed 4x) to identify pure components. Model was clean and easy to obtain. Then CLS was performed on each individual image at full spatial resolution (RAPID) to obtain final concentrations of each component.

MCR Analysis




This is the percent variance averaged across all of the data – this varies from image to image, depending on the image.

Questions: Spectrally I would interpret Factor 4 to be DAPI (we collect only the tail of it), however if Dapi stains the cell nuclei this is not correct. Neither factor 3 or 4 really look right for DAPI .

Summary


GFP




Min = 0
Max = 1000

If this were on the same
scale as the others this
would be black, model is
just fitting noise.

Control, GFPexp




CCL2, GFPexp





Min = 0
Max = 31000

AF red




Min = 200
Max = 16000





Min = 200
Max = 6000

AF teal



Min = 500
Max = 12000




Min = 500
Max = 2500

C5781_6_control_1_4.dat

CX3_CR1_GFP_Control_1_9.dat


CX3_CR1_GFP_INJ_CCL2_1_9.dat

Sample RGB Images




Control/No GFP

C5781_6_control_1_4.dat



Control/GFP exp

CX3_CR1_GFP_Control_1_9.dat




CCL2/GFP exp

CX3_CR1_GFP_INJ_CCL2_1_9.dat


Note: Colorscale is NOT the same for each image.

Sample RGB Images




Control/No GFP

C5781_6_control_1_3.dat



Control/GFP exp

CX3_CR1_GFP_Control_1_5.dat



CCL2/GFP exp

CX3_CR1_GFP_INJ_CCL2_1_5.dat

Note: Colorscale is NOT the same for each image.



These next slides show the individual components for all images of CC2 treated and untreated retinas

All images on the following slides are 100um x 100um and all images within a slide are plotted on the same scale to facilitate direct comparison and assessment of morphological differences.



Note: I don't know why the bottom right image is messed up. The original looks fine. My guess is something happened in the preprocessing since we are right at the edge of the memory available.

GFP expressing - Control

This image was
run at OD 2.5,
rather than 1 like
pthers.




GFP- Expressing – CCL2




Component 2 = GFP

Control – GFP expressing

This image was
run at OD 2.5,
rather than 1 like
pthers.






GFP- Expressing – CCL2





Control – GFP expressing

This image was
run at OD 2.5,
rather than 1 like
pthers.



GFP- Expressing – CCL2



Component 4 AF (diffuse)

Summary

- What are the identities of the two AF components? Is one a dye?
- No drastic differences (intensity, morphology, or spectral shift) between control and CCL2 treated retinas in any component (GFP or AF) in the GFP expressing mice. Subtle differences are possible...
 - Maybe a slight hint of the processes being retracted in the CCL2 treated images
 - AF (blue) is possibly more intense in the CCL2 treated mice
- There's an intensity difference in the non-GFP data AF components compared to the GFP expressing AF components, however this is due to different excitation intensities during the scan.