


### We extract the surface potential from photoemission spectra





## We need ionization energy to associate the surface potential to the electronic properties of semiconductors

Threshold photoemission & ionization energy









#### Surface potential & ionization energy match well in metallic materials

Example of a metallic sample

10 μm

4000

3000

2000

1000

-1.0

-0.8

-0.6

Relative work function (eV)

Number of pixels

- Multi-layer graphene on a SiC substrate displays varying work function depending on the layer thickness
- The distribution width is ~50meV consistent to kT for RT

5000 -


4000 -

3000 -

2000 -


1000


4.0





### The surface potential & ionization energy vary in the vicinity of v-defects

Defects' locations are identified unequivocally in LEEM image









Small v-defects & trench defects have ~40meV
lower surface potential than the pristine area

- 0.0


- -0.1


eV

Large v-defects are ~150meV lower


## Our surface potential & ionization energy measurements implies the energy barrier for holes at v-defect site

First, define the spatial variation of the surface potential ( $E_F$  is unknown)






Determine the energy of the VBM w/ respect to the surface potential




 Significant barrier heights (possibly > 0.3 eV) for holes can be expected at the vicinity of v-defects & trench defects

Define the energy of CBM based on VBM for the pristine area

Estimate the electronic states for InGaN MQW



## Our experiment suggests defect screening in MQW without compositional variation at sub-micron length-scale




Potential

"Defect Screening

carriers in QWs

Lateral Position

# Perspectives of LEEM-PEEM research and instrumentation in the bigger semiconductor research



- Laser-based PEEM instrumentation for higher intensity and tunable photon source
  - Addressing issues related to microscope's alignment and stability
  - Use of various optical transitions to access various energy levels of semiconductors
  - Time- & spin-resolution to address the fundamental questions for materials envisioned for quantum computing
- In-operando microscopy to track device operation, chemical reaction, etc.



Sandia

"Sandia National Laboratories is a multi-program laboratory managed and operated by Sandia Corporation, a wholly owned subsidiary of Lockheed Martin Corporation, for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-AC04-94AL85000."